

voué

# Platinum NIR Milk Light







# Rejuve

Skin Rejuvenation . Rejuve your beauty


# 7 MAJOR EFFECTS

Rejuvenate the skin

## The first near-infrared light wave technology provider in China.

The NIR Milk Light uses a high power light source in the near infrared spectrum to achieve safe and effective dermal heating, imparting a youthful glow to the skin. The pleasant treatment elevates dermal temperature, generates thermal injury to the tissue which contracts existing fibers, stimulates the formation of new collagen and improves its alignment and thickness

"Milk light", that is, near-infrared light (NIR) is between visible light and infrared light, using near-infrared technology, the selected spectrum range is 900-1600nm, and the peak is concentrated in the infrared light wave of 1300nm. The target tissue is the water molecule in the skin. NIR infrared light waves can penetrate the epidermis and reach the dermis layer, effectively heating the water molecules in the dermis, promoting the contraction of collagen fibers and stimulating the production of new collagen fibers. Its penetration depth in the skin is the deepest of all light and lasers, much deeper than photons.



Improve skin tone

1

Shrink pores

2

Fight anti-aging

3

Lighten fine lines

4

Lighten stains

5


Reduce redness

6


Stimulate collagen  
regeneration

7





# Platinum NIR Milk Light



**Innovation:** Platinum NIR near-infrared light & ssr narrow-spectrum light the world's unique precision light technology.

**Activation:** REJUVE NIR 1300nm wave near-infrared light, really empowers cells to promote deep collagen repair and regeneration.

**Depth:** Effectively reach the skin 0.5 ~ 2mm, directly to the middle layer of the dermis to the deeper layer of the skin.

**Function:** Revitalizes the basement skin layer with connective tissue regeneration, and restores healthy skin.



## Doing this NIR Milk Light,

*Embrace Youthful Glow to the Skin*

**Skin texture + Complexion full effect management.**

The only narrow-spectrum light with a wavelength of only 100nm in the world.

## Unique sliding technology

Ultra-high frequency + fast sliding technology brings a painless, comprehensive, safe and comfortable experience, and makes the skin fair and translucent immediately.

## Heraeus particular kind infrared light source

The light source comes from the special light source manufacturer - Heraeus Company of Germany, with 1000 hours of energy to ensure stable work and safer for the skin.







### Therapeutic advantage

Non-invasive, painless, no recovery period.  
No scab, redness, inflammation and other side effects.  
Instantly whiten and brighten the skin instantly.  
The treatment process is efficient and comfortable.

### Application advantages



#### Precise light wave

effective depth of 0.5 ~ 2mm, rejuvenate cells.



#### Double-Handle and double spot

face + body one-stop full-body rejuvenation.



#### Comfortable treatment

stronger and more comfortable than traditional photon effects.







#### Cooling technology

Intelligent cooling system, humanized skin protection mechanism, bring better experience.

### Application

People with skin problems such as dry and rough skin, pigmentation, skin sensibility, dull complexion, acne, pimples, pox marks, photoaging, cutis laxa, etc.



### Why to choose Platinum NIR Milk Light?

Originated from Israel, hundreds of scientists and doctors have 15 years of technology research and development and clinical practice to create Platinum NIR.

1300nm wave peak near-infrared light, continuously revitalize the basement membrane belt, restore healthy skin, accelerates melanin metabolism and brightens skin tone.



Before

After



Before

After




Before



Before



After



After